

Instruction Manual for the Ray *de novo* genome assembler software

Sébastien Boisvert

February 23, 2011

Ray version 1.2.4

<http://denovoassembler.sf.net>

Reference to cite:

Sébastien Boisvert, François Laviolette & Jacques Corbeil.

Ray: simultaneous assembly of reads from a mix of high-throughput sequencing technologies.

Journal of Computational Biology (Mary Ann Liebert, Inc. publishers, New York, U.S.A.).

November 2010, Volume 17, Issue 11, Pages 1519-1533.

doi:10.1089/cmb.2009.0238

<http://dx.doi.org/doi:10.1089/cmb.2009.0238>

Contents

1	Installation	3
1.1	Compilation flags	3
1.2	Compilation with the configure script	3
1.3	Compilation with the makefile.alternate	3
2	Inputs	3
3	Parameters	3
3.1	k -mer length with -k	3
3.2	Output prefix with -o	4
3.3	Single-end reads with -s	4
3.4	Paired-end reads and mate-pair reads with -p	4
3.5	Paired-end reads and mate-pair reads with -i (interleaved sequences)	4
4	Output	4
4.1	Contiguous sequences	4
4.2	Paired-end and mate-pair libraries	4
4.3	Coverage distribution	5
4.4	Messages	5
5	Example	5
5.1	Bacterial genome with paired-end and mate-pair short reads . . .	5

1 Installation

1.1 Compilation flags

For processors requiring addresses aligned on 8 bytes, you must add the following compilation flag.

```
-DALIGN_ADDRESSES
```

1.2 Compilation with the configure script

```
./configure --prefix=/software/ray-x.y.z
make
sudo make install
ls -l /software/ray-x.y.z/bin/Ray
```

1.3 Compilation with the makefile.alternate

```
make -f makefile.alternate
ls -l code/Ray
```

2 Inputs

The input files for Ray contain sequences. The files must be formatted in one of the supported formats.

Table 1: File formats compatible with the Ray *de novo* genome assembler software

Format	Obligatory extension
Fasta format	.fasta
Fasta format, compressed with GNU zip (gzip)	.fasta.gz
Fasta format, compressed with bzip2	.fasta.bz2
Fastq format	.fastq
Fastq format, compressed with GNU zip (gzip)	.fastq.gz
Fastq format, compressed with bzip2	.fastq.bz2
Standard flowgram format	.sff

3 Parameters

3.1 k -mer length with -k

Ray builds a distributed catalog of all occurring k -mers in the reads and their reverse-complement. k must be greater or equal to 15 and lower or equal to 32.

3.2 Output prefix with -o

Output files are named according to the prefix provided by the option -o.

3.3 Single-end reads with -s

-s <sequencesFile>

3.4 Paired-end reads and mate-pair reads with -p

-p <leftSequencesFile> <rightSequencesFile> [<averageFragmentLength> <standardDeviation>]

Example for paired-end reads (the ends of DNA fragments):

-p s_200_1.fastq s_200_2.fastq

Example for mate-pair reads (owing to circularization and ligation, the order must be reverse):

-p s_20000_2.fastq s_20000_1.fastq

3.5 Paired-end reads and mate-pair reads with -i (interleaved sequences)

-i <sequencesFile> [<averageFragmentLength> <standardDeviation>]

In the interleaved file (example is for a fasta file):

```
>200_1_1234/1
ATCGATCGATCGACTCAGACACGTACG
>200_1_1234/2
ACTGACGACGTACGACGTCATGCAACT
...
```

For mate-pair reads, the order must be reverse.

4 Output

4.1 Contiguous sequences

OutputPrefix.fasta contains contiguous sequences.

4.2 Paired-end and mate-pair libraries

OutputPrefix.LibraryLibraryNumber.txt contains the distribution of distances for paired-end and mate-pair libraries. One file per library.

4.3 Coverage distribution

OuputPrefix.CoverageDistribution.txt contains the k -mer coverage distribution.

4.4 Messages

OuputPrefix.ReceivedMessages.txt contains a matrix.

It contains the number of received messages for each MPI rank.

(MPI communication matrix; rows=destinations, columns=sources)

5 Example

5.1 Bacterial genome with paired-end and mate-pair short reads

The command:

```
mpirun -np 32 ~/Ray/trunk/code/Ray \  
-p /home/boiseb01/nucore/Large-Ecoli/200_1.fastq \  
  /home/boiseb01/nucore/Large-Ecoli/200_2.fastq \  
-p /home/boiseb01/nucore/Large-Ecoli/1000_1.fastq \  
  /home/boiseb01/nucore/Large-Ecoli/1000_2.fastq \  
-p /home/boiseb01/nucore/Large-Ecoli/4000_1.fastq \  
  /home/boiseb01/nucore/Large-Ecoli/4000_2.fastq \  
-p /home/boiseb01/nucore/Large-Ecoli/10000_1.fastq \  
  /home/boiseb01/nucore/Large-Ecoli/10000_2.fastq \  
-o BacterialGenome | tee RayLog
```